

Glass Syringes

TRÜTH Interchangeable Glass Syringes

Ideal for laboratory liquid handling of various liquids including Plastic Sensitive Materials, hazardous chemicals, Adjuvant solutions, petroleum derivatives, glues, synthetic dyes, perfumes, essential oils and organic chemicals.

- Glass Barrels & Pistons are made from borosilicate glass provide protection against thermal shock
- Graduation marks are permanently fused into a glass barrel
- Standard Luer Lock Syringes confirm to ISO 594/1 ; 595/2
- Precision shrinking ensures improved leak free smooth operation
- Coloured plunger facilitates dosage measurements
- Autoclavable at 121°C
- We supply in private label terms as well
- Luer locks are made of Brass Nickel-Chrome & Stainless Steel (SS316) plated
- Inert liquid path, No 'O' ring



Termination	Volume	Graduation	Packing	Catalogue
Center all	1ml	0.05ml	10pcs	01-03-01-03
Glass Tip	2ml	0.1ml	10pcs	01-05-01-04
(CT)	3ml	0.1ml	10pcs	01-07-01-04
	5ml	0.2ml	10pcs	01-08-01-05
	10ml	0.2ml	10pcs	01-09-01-05
	20ml	0.5ml	10pcs	01-10-01-06
	30ml	1ml	1pc	01-11-01-07
	50ml	2ml	1pc	01-13-01-08
	100ml	2ml	1pc	01-15-01-08

Termination	Volume	Graduation	Packing	Catalogue
Metal	1ml	0.05ml	10pcs	01-03-02-03
Luer	2ml	0.1ml	10pcs	01-05-02-04
Lock Tip	3ml	0.1ml	10pcs	01-07-02-04
(CL)	5ml	0.2ml	10pcs	01-08-02-05
	10ml	0.2ml	10pcs	01-09-02-05
	20ml	0.5ml	10pcs	01-10-02-06
	30ml	1ml	1pc	01-11-02-07
	50ml	2ml	1pc	01-13-02-08
	100ml	2ml	1pc	01-15-02-08

Termination	Volume	Graduation	Packing	Catalogue
Stainless	1ml	0.05ml	10pcs	01-03-75-03
Steel (SS316)	2ml	0.1ml	10pcs	01-05-75-04
Luer Lock Tip	3ml	0.1ml	10pcs	01-07-75-04
(SS-CL)	5ml	0.2ml	10pcs	01-08-75-05
	10ml	0.2ml	10pcs	01-09-75-05
	20ml	0.5ml	10pcs	01-10-75-06
	30ml	1ml	1pc	01-11-75-07
	50ml	2ml	1pc	01-13-75-08
	100ml	2ml	1pc	01-15-75-08

Also available with Metal Tip  on request