

Storage Vials for Custom Synthesis & Impurities

1) 0.5ml Micro Storage Vial, (U-Tube Vial)

Catalogue No : TC1650A

Description : Amber Vial with Blue Closed Cap & liner

German  Technology



MATERIALS

Vial:
Borosilicate
Glass, Type 1


Screw Cap:
Polypropylene

Septa / liner:
Silicone with
PTFE

Other Storage Vials with Black Cap + Septa / Pack of 100

Sr.No	Size	Catalogue No	OD x Height
2	2ml, Clear	TC-80A + TC-804A	11.6 x 32mm
3	, Amber	TC-80C + TC-804A	
4	4ml, Clear	TC-130A + TC-1304K	14.7 x 45mm
5	, Amber	TC-130C + TC-1304K	
6	20ml, Clear	TC-2420A + TC-2404L	27.5 x 57mm
7	, Amber	TC-2420C + TC-2404L	
8	30ml, Clear	TC-2430A + TC-2404L	27.5 x 72.5mm
9	, Amber	TC-2430C+ TC-2404L	
10	40ml, Clear	TC-2440A + TC-2404L	27.5 x 95mm
11	, Amber	TC-2440C+ TC-2404L	
12	60ml, Clear	TC-2460A + TC-2404L	27.5 x 140mm
13	, Amber	TC-2460C+ TC-2404L	

5ml & 10ml Storage Vial (Coming Soon)
 Amber Vial with Plastic Cap + PP Screw Cap



ADDRESS:

Top Syringe Compound, W.E. Highway, Pandurang Wadi,
 Mira Road (East), Thane 401 104. INDIA
 Tel.+91 22 6520 8849

Email: sales@top-syringe.com

Website: www.top-syringe.com